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## **NEMA Standards Publication NU 2-2018**

*Performance Measurements of Positron Emission Tomographs (PETS)*

*Itself based on*

## **IEC 61675-1:2013**

Radionuclide imaging devices -  
Characteristics and test conditions - Part 1:  
Positron emission tomographs

## **NEMA Standards Publication NU 1-2018**

*Performance Measurements of Gamma Cameras*

*Itself based on*

## **IEC TR 61948-2:2001**

Nuclear medicine instrumentation - Routine tests - Part 2: Scintillation cameras and single photon emission computed tomography imaging


These Standards Publications establish definitions, quantitative measurements of performance characteristics, and reporting techniques for the specification of the various camera parameters.

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The European Association of Nuclear Medicine (EANM), with its central role in the scientific community and its professional setup, today constitutes the umbrella organisation for nuclear medicine in Europe. In 2006, the leadership of the EANM launched EANM Forschungs GmbH (EARL) as an initiative to promote multicentre nuclear medicine and research.

## ACR–AAPM TECHNICAL STANDARD FOR MEDICAL NUCLEAR PHYSICS PERFORMANCE MONITORING OF PET IMAGING EQUIPMENT

### III. PERFORMANCE CHARACTERISTICS TO BE MONITORED

#### A. Acceptance Testing

Initial performance testing of imaging equipment must be performed upon installation and should be completed before clinical use. This testing should be more comprehensive than periodic performance testing and should be consistent with current acceptance testing practices [4-8].

#### B. Performance Evaluation


The Qualified Medical Physicist must design a quality assurance (QA) program that includes regular testing procedures to insure proper operation on a daily basis. The PET testing procedures must follow the manufacturer's prescribed schedule and must include, at a minimum, daily quality control (QC) tests as well as quarterly scanner testing to evaluate scanner calibration, uniformity, resolution, and contrast. The QA program must be reviewed at least annually by the Qualified Medical Physicist who must be onsite at least annually.

The following characteristics should be evaluated on at least an annual basis as applicable to the design of the scanner [8,9]:

1. Spatial resolution
2. Count rate performance (count rate versus activity), including count loss correction
3. Sensitivity (cps/MBq/mL)
4. Image uniformity
5. Image quality
6. Accuracy of attenuation and scatter correction, and SUV measurement
7. Safety evaluation
  - a. Mechanical
  - b. Electrical

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






Harmonization based on quantitative analysis of standardized uptake value variations across PET/CT scanners: a multicenter phantom study  
Monsef, Abbasa,b; Ay, Mohammad Rezaa,b; Sheikhzadeh, Peymana,c; Geramifar, Parhamd; Rahmim, Armane,f; Ghafarian, Pardisg,h; Nuclear Medicine Communications 43(9):p 1004-1014, September 2022. | DOI: 10.1097/MNM.0000000000001598

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Harmonization based on quantitative analysis of standardized uptake value variations across PET/CT scanners: a multicenter phantom study  
 Monsef, Abbasa,b; Ay, Mohammad Rezaa,b; Sheikhzadeh, Peymana,c; Geramifar, Parhamd; Rahmim, Armane,f; Ghafarian, Pardisg,h; Nuclear Medicine Communications 43(9):p 1004-1014, September 2022. | DOI: 10.1097/MNM.0000000000001598

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
Peters et al., *Towards standardised multi-centre and multi-vendor SPECT/CT* [9]

Coming up soon  
**<sup>177</sup>Lu SPECT/CT**  
Accreditation

**Lu-177 EARL Quantitative SPECT**  
prototype  
15Mar2022

©2022 Research4Life, Inc.  
All rights reserved. No part of this document may be reproduced without prior written permission from Research4Life, Inc.  
19100 No. W. Avenue, #100, Denver, CO 80233

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



FIG. 27. Image of a system spatial resolution phantom consisting of two sets of parallel line sources filled with  $^{99m}\text{Tc}$  in solution. The solution is contained in glass microcapillary tubes, about 0.5 mm in internal diameter, plugged at each end with sealing clay or wax and mounted on the edges of a 10 cm square drawn on a sheet of rigid plastic or heavy cardboard. Note that the resolution in the X and Y directions differ.

6 Concentric aluminium tubes used to detect camera sensitivity in PET  
Performance Measurements of Scintillation Cameras, NEMA Standards  
Publication No. NU2, National Electrical Manufacturers Association (NEMA),  
Washington, D.C., 2001

**IAEA NMQC  
Toolkit**

Free ImageJ plugins to support nuclear  
medicine physics equipment quality control  
examinations.

<https://humanhealth.iaea.org/HHW/MedicalPhysics/NuclearMedicine/QualityAssurance/NMQC-Plugins/index.html>

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The **user** is responsible for performing and documenting the stability checks. In particular, in order to ensure the independence of these measurements, the user **cannot mandate the supplier** or the installing company to perform the stability checks. If a stability check fails, the user has to contact the responsible **medical physicist** and corrections must be made before the machine can be used clinically.

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# Dokumentation

Content	Remarques
Indication of the FOPH licence number	The FOPH license number for the nuclear medicine camera concerned must be indicated (for SPECT-CT and PET-CT devices: number of the license to operate a radiology system for medical use; for independent gamma cameras (without CT): number of the license to use unsealed radioactive sources).
Reception test reports	Includes records of reception tests and records of handover of the installation to the user. Target values and tolerance intervals, if not specified in the Annexes to this guideline, must correspond to the manufacturer's specifications and must be documented and passed on to users.
Status checks reports	Includes records of status checks. Target values and tolerance intervals, if not specified in the Annexes to this guideline, must correspond to the manufacturer's specifications and must be documented and passed on to users.
Revision checklists (maintenance) and periodicity	According to the manufacturer's data.
<b>Procedures and reference values for the stability checks</b>	Certain <b>reference values for later stability checks are defined at the time of the reception test</b> or are redefined at the time of each status check. These reference values must be documented. The procedures to be followed for the stability checks must also be defined by the manufacturer. The reports of the stability checks are filed separately, either in electronic form or on paper.
Indication of the version and storage location of the instructions for use	The instructions for use must be available at all times in the language used in the unit.
CE Declaration of Conformity	In accordance with the Medical Devices Ordinance (MedDO; SR 812.213; [10]) a CE Declaration of Conformity must be available for the device.
Indication of ambient dose rates in the room where the camera is installed and in adjacent rooms	The maximum permitted ambient dose rates must be respected in the camera room and in all adjacent rooms. For PET equipment, the radiation protection plans and calculation tables in accordance with the guideline "Radiation protection measures relating to the construction of PET installations" [11] must be available.

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N°	Parameter	Periodicity	Reference and method	Devices concerned	Nuclids, activities and phantoms	Associated stability check
G1	Intrinsic flood field uniformity	RT and SC	NEMA NU 1-2018, part 2.4 The intrinsic uniformity of the system is measured, i.e., the response of the system without a collimator to a uniform flux of radiation from a point source.	All, except those ones for which the collimators cannot be removed	Contrarily to what is requested by NEMA, it is sufficient to measure this parameter for two nuclides, Tc-99m and another one with photons of higher energy (except if only Tc-99m is used).	--
G2	System image homogeneity & collimator check	RT and SC	Both the condition of the collimators and the images produced with the collimators (homogeneity) shall be checked visually.	All	RT: check of all the collimators. SC: check of the mostly used collimator.	--
G3	Intrinsic energy resolution	RT and SC	NEMA NU 1-2018, part 2.3 The ratio of the photopeak FWHM to the photopeak center energy shall be calculated.	All	Use of Tc-99m and/or Co-57 if the latter is used for daily quality assurance. Note that Tl-201 shall not be used as a reference isotope.	--
G4	Intrinsic spatial resolution	RT and, if KG5 is out of tolerance, SC	NEMA NU 1-2018, part 2.1 The FWHM and FWTM shall be determined.	All, except CZT detectors (Cadmium-Zinc-Telluride)	Use of Tc-99m.	--
G5	System spatial resolution without scatter	RT	NEMA NU 1-2018, part 3.1 The FWHM and FWTM of the line spread function shall be measured.	All	The test has to be performed for all collimators. Use of two capillary tubes.	KG5
		RT SC	Visual check of the spatial resolution. If the quadrant bar phantom is used for KG5, a reference value (best resolution visible) shall be determined during RT.	All	The test has to be performed for all collimators. Use of a dedicated quadrant bar phantom.	KG5
G6	System planar sensitivity and collimator penetration and scatter	RT and SC	NEMA NU 1-2018, part 3.3 The ratio of collimated counts to the known activity of a planar source shall be determined.	All	RT: all the collimators. SC: with the mostly used one only. The activity has to be injected into Petri dishes (as requested by NEMA), but the use of traceable sources is allowed.	--
			NEMA NU 1-2018, part 4.4 Measurement of the total system sensitivity to a uniform concentration of		Use of Tc-99m and the nuclides used for dosimetry (which includes planification of	

			visible) shall be determined during RT.			
G6	System planar sensitivity and collimator penetration and scatter	RT and SC	NEMA NU 1-2018, part 3.3  The ratio of collimated counts to the known activity of a planar source shall be determined.	All	RT: all the collimators. SC: with the mostly used one only.  The activity has to be injected into Petri dishes (as requested by NEMA), <b>but the use of traceable sources is allowed.</b>	--
	System volume sensitivity	SC RT	NEMA NU 1-2018, part 4.4  Measurement of the total system sensitivity to a uniform concentration of activity, and determination of the average volume sensitivity per axial centimeter.  For those devices that need to perform the corresponding stability checks KG4, reference values shall be determined during RT.	SPECT devices	Use of Tc-99m and the nuclides used for dosimetry (which includes planification of the activity to be administrated for a therapy).  Use of a specific NEMA fillable cylindrical phantom.	KG4
G7	System count rate performance with scatter	RT	NEMA NU 1-2018, part 3.5 or AAPM report 117  The maximum observed count rate and the observed count rate for a 20% count loss shall be measured, and a curve of observed count rate with scatter versus input count rate shall be made.  As an alternative to the NEMA methodology, the use of one of the methods mentioned in the AAPM report 117 [12] is allowed.	All	Use of Tc-99m and of a cylindrical phantom.	--
G8 Vormals G9	Intrinsic spatial linearity	RT	NEMA NU 1-2018, part 2.2  The goal of the test is to identify potential distortions in the images. The intrinsic spatial differential and absolute linearity has to be reported.	All, except CZT detectors	Use of Tc-99m and a lead mask with parallel slits.	--
	System linearity	RT <del>SC</del> -> KP	In addition, a visual check for linearity and absence of distortions shall be performed.	All	Use of the data from G5 (both parts).	KG6
G9 Vormals G11	Whole-body system spatial resolution without scatter	RT <del>SC</del> -> KP	NEMA NU 1-2018, part 5.1  The system spatial resolution without scatter shall be measured parallel and perpendicular to the direction of continuous motion and expressed as FWHM or FWTM of the line spread function.	All, except step-and-shoot whole-body planar acquisition devices	Use of two capillary tubes of Tc-99m. The camera shall be equipped with a collimator.	--
			NEMA NU 1-2018, part 4.1, or according to the manufacturer			

<p><b>G10</b> <b>Vormals G12</b></p>	<p>System alignment</p>	<p>RT and SC</p>	<p>or the line spread function. NEMA NU 1-2018, part 4.1, or according to the manufacturer</p> <p>The transaxial alignment of acquired images with the system's mechanical center of rotation shall be measured. Likewise, for multi-head SPECT imaging systems, the axial alignment of images from the individual heads shall be measured.</p>	<p>SPECT devices only</p>	<p>If the NEMA method is used, three Tc-99m or Co-57 point sources shall be used.</p>	<p>KG7</p>
<p><b>G11</b></p>	<p>Detector-detector sensitivity variation</p>	<p>RT</p>	<p>NEMA NU 1-2018, part 4.5</p> <p>Assessment of the relative difference in sensitivity of the individual detectors. All the projection images collected using each detector shall be summed up. The difference between the largest and the smallest of these quantities shall be computed.</p> <p>Reference values shall be determined for the later KG9.</p>	<p>All multi-detector camera systems, when applicable</p>	<p>Use of the data from the second part of G6. The use of one collimator is enough.</p>	<p>KG9</p>
<p><b>G12</b></p>	<p>Tomographic contrast and absolute quantification accuracy</p>	<p>RT</p>	<p>NEMA NU 1-2018, part 6</p> <p>Images simulating those obtained in a total body imaging study with hot and cold lesions of different diameters are produced. Measurement of the contrast of cold and hot spheres in a warm background, of the variability of the background regions, and of the deviation of the large, cold reconstructed lung region without activity.</p> <p>For centers using absolute quantification (image output in Bq/ml), the absolute activity concentration shall also be measured and compared to the injected activity concentration.</p> <p>Reference values shall be determined for the later KG8.</p>	<p>SPECT devices only</p>	<p>Use of the NEMA/IEC image quality phantom (same phantom as for the PET reception test P5). Use of Tc-99m. Background activity shall be 20 kBq/ml. The sphere-to-background activity concentration shall be 8:1.</p>	<p>KG8</p>
<p><b>G13</b></p>	<p>Image homogeneity and precision of quantification</p>	<p>RT</p>	<p>Verification of the background activity. The mean coefficient of variation (image homogeneity) and mean background activity concentration (precision of quantification) shall be computed according to the paragraph "Analysis" of</p>	<p>Image homogeneity: all devices.</p> <p>Precision of quantification: only centers undertaking quantitative studies or dosimetry</p>	<p>Use of the data from the second part of G6.</p> <p>The test has to be performed for the same nuclides as those defined in the second part of G6. Tolerance interval: see paragraph</p>	<p>KG10</p>

G13	Image homogeneity and precision of quantification	RT	<p>the later KG8.</p> <p>Verification of the background activity. The mean coefficient of variation (image homogeneity) and mean background activity concentration (precision of quantification) shall be computed according to the paragraph "Analysis" of Annex 6.</p>	<p>Image homogeneity: all devices.</p> <p>Precision of quantification: only centers undertaking quantitative studies or dosimetry.</p>	<p>Use of the data from the second part of G6.</p> <p>The test has to be performed for the same nuclides as those defined in the second part of G6. Tolerance interval: see paragraph "Tolerance" of Annex 6.</p>	KG10
G14	SPECT-CT coregistration accuracy	RT and SC	<p>Preferably NEMA NU-1 2018, part 7, or according to the manufacturer</p> <p>NEMA NU-1 2018 method: The alignment accuracy between SPECT and CT image volumes shall be measured, by using data acquired with SPECT and CT fiducial markers at six locations.</p>	SPECT-CT only	<p>Radionuclide for the SPECT portion of the fiducial markers: Tc-99m.</p> <p>It is allowed not to use the weights prescribed by NEMA NU-1 2018.</p>	--

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

## System Spatial Resolution

Quadrant Bar Phantom or Capillary Phantom:

N°	Parameter	Periodicity	Reference and method	Devices concerned	Nuclids, activities and phantom
KG1	Background count rate	At the beginning of each working day	Perform a blank scan of 2.5 to 3 minutes. A given number of counts (manufacturer-specific) should not be exceeded to make sure that no radioactive contamination is present in the system. In case of contamination, measures have to be taken.	All	--
KG2	Check of the energy window	Qualitatively, before each exam Quantitatively, once a week	Check of the correct position of the Tc-99m peak. A visual check shall be performed before each exam, with the injected patient himself.  Once a week, a quantitative measure shall be performed and documented, in order to follow any shift in the energy window.	All	Use of Tc-99m.
KG3	Homogeneity	Half-yearly Daily !!!	Methodology (with or without collimator, maximal number of counts, point source or flat source) according to the manufacturer.	All	See column "Reference and method"
KG4	System volume sensitivity	Yearly	The same method as described in the second part of G6 shall be followed and the results shall be compared to the reference values of the reception test.	Quantitative SPECT devices	See G6. The data collected for this test used as well for KG9 and KG10.
KG5	System spatial resolution without scatter	Yearly	Use of one of the two methods of G5. If the quadrant bar phantom is used, the results shall be compared to the reference values of the G5 reception test.	All	The test has to be performed with collimators.  The test shall be performed either with the quadrant bar phantom for measuring spatial resolution (visual check) or with capillary tubes (determination of FWHM and FWTM values).
KG6	System linearity	Yearly	Visual check of the absence of distortions in the images.	All	Use of the same lead mask as in the first part of G8, or the quadrant bar phantom used for measuring spatial resolution, and Tc-99m or a Co-60 source. Alternatively, Tc-99m capillary tubes can be used; see [11]. The same

# SPECT: KoI

KG6 System linearity

## Annex 6 Detailed instructions for reception tests G13 and P7 and stability checks KG10 and KP6: image homogeneity and precision of quantification

**Purpose:** To verify the accuracy of image quantification and image homogeneity.

**Materials:** **Cylindrical phantom** which should cover the axial FOV filled with 18F. Diameter of the phantom should be at least 20 cm.

KG7 System alignment

**Procedure:** Prepare a syringe with about 50 - 100 MBq of 18F by accurately measuring it in the same dose calibrator used for clinical routine, writing value of activity and time of assay. The clocks used for recording the time assay should be checked against the scanner time. Put the activity in the phantom, mixing with water solution thoroughly to get a uniform radioactivity distribution. Measure the residual activity in the syringe and evaluate the net activity put in the phantom. Add water until the phantom is completely filled. Place the phantom on the phantom holder provided by the manufacturer and move the phantom at the beginning of the CT FOV, making sure that the phantom is centered in the FOV also with respect to the height. In case of no phantom holder, simply place it on the patient's pallet.

KG8 Tomographic contrast and absolute quantification accuracy

KG9 Detector-detector sensitivity variation

The acquisition should be performed by using the protocol provided by the manufacturer (**same protocol for the reception tests P7 and the later stability checks KP6**). In case no acquisition protocol is provided, the standard protocol used for clinical routine should be used (body, head, ...). Make sure that a low dose CT for attenuation and scatter correction purposes is included in the procedure. Attention should be paid when inserting data regarding 18F activity assayed by the dose calibrator, time of assay and weight of the volume of the 18F solution used to fill the phantom. Acquisition stop condition should be set for 100 million counts.

KG10 Image homogeneity and precision of quantification

**Analysis:** For the acquired images, draw 1 circular ROI with a diameter greater than 5 cm, at least 2 cm from the phantom edge on the reconstructed central slice and on  $\pm 5$ -6 adjacent slices (total of 12 ROIs). Measure the mean and standard deviation activity concentration for each ROI and calculate the coefficient of variation (COV) for each slice:

$$\text{coefficient of variation (\%)} = \frac{\text{standard deviation}}{\text{mean}} \times 100$$

The mean COV and mean background activity concentration can then be calculated among the different ROIs.

**Tolerance:** The ideal mean background activity concentration should be within  $\pm 5\%$  of the true activity concentration and the mean coefficient of variation should be  $< 15\%$ .

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N°	Parameter	Periodicity	Reference and method	Devices concerned	Nuclids, activities and phantoms	Associated stability check
P1	Spatial resolution	RT	NEMA NU 2-2018, part 3 Imaging point sources in air and then reconstructing images with no smoothing or apodization . The FWHM and FWTM of the image point spread function are measured.	All	Use of F-18 or Na-22.	--
P2	Scatter fraction, count losses, and randoms	RT	NEMA NU 2-2018, part 4 The relative system sensitivity to scattered radiation, the effects of system dead time and the generation of random events are measured.	All	Use of a dedicated phantom (solid right circular cylinder with a hole drilled parallel to the central axis of the cylinder) with approximately 500 MBq of F-18.	--
P3	Sensitivity	RT	NEMA NU 2-2018, part 5 The sensitivity of the system (rate in counts per second that true coincidence events are detected for a given source strength) is measured.	All	Use of F-18 with approximately 10 MBq activity. Use of the NEMA sensitivity phantom (five sleeves of different diameters). The use of a cylinder with Ge-68 sources (according to the former NEMA NU 2-1994) is not allowed.	--
P4	Accuracy: corrections for count losses and randoms	RT	NEMA NU 2-2018, part 6 The accuracy of corrections for dead time losses and random event counts is measured.	All	Test performed with the measurement data from P2.	--
P5	Image quality, accuracy of corrections AC-Correction here AC-Map (former P8)	RT	NEMA NU 2-2018, part 7 Images simulating those obtained in a total body imaging study with hot lesions of different diameters are produced, with activity also present outside the region of interest of the scanner to reproduce the clinical routine. Image contrast and background variability ratios for hot spheres are used as measures of image quality. The accuracy of corrections is determined from the uniform background and cold lung insert regions.	All	Use of the NEMA/IEC image quality phantom (body phantom). Radionuclide to be used: F-18 with 370 MBq for 70 kg. Use as well of the solid right circular cylinder (scatter phantom, same as for P2) with same activity concentration. Fill the spheres with a unique activity ratio of 4:1.	--
	Qualitative contrast		Qualitative determination of the contrast by counting and reporting the number of visible lesions in images simulating those obtained in a total body imaging study with hot lesions of different diameters		Use of the NEMA/IEC image quality phantom prepared for P5 one hour after P5 has been performed. The test has to be performed <u>without</u> the scatter phantom and with the clinically recommended algorithm.	

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## Quantitative contrast test

NEMA/IEC Image quality phantom

N°	Parameter	Periodicity	Reference and method	Devices concerned	Nuclids, activities and phantom	
KP1	PM check	Every working day	Check of the amplification factor (gain), offset PM, and homogeneity.	All	Use of an external positron source, except devices equipped with an appropriate internal source (for example lutetium sources).	
KP2	Control of the energy window	Every working day	Check of the setting and FWHM resolution.	All	Use of an external positron source, except devices equipped with an appropriate internal source (for example lutetium sources).	
KP3	Coincidence timing	Weekly check for the semiconductor devices, on every working day for the other (PM) devices.	The parameter to be checked is as defined by the manufacturer.	All	Use of an external positron source.	
KP4	Visual verification of the system	Every working day	Comparison of the sinograms.	All	Use of an external positron source, except devices equipped with an appropriate internal source (for example lutetium sources).	According to the manufacturer <sup>12</sup>
KP5	Quantitative test Qualitative contrast test	Yearly	<b>Qualitative determination of the contrast</b> by counting and reporting the number of visible lesions in images simulating those obtained in a total body imaging study with hot lesions of different diameters.	All	Prepare the NEMA/IEC image quality phantom with 250 MBq of F-18 and fill the spheres with a unique activity ratio of 4:1. The test has to be performed <u>without</u> the scatter phantom and with the clinically recommended algorithm.  For the centers undergoing accreditation tests for quantitative measurements (e.g., EARL), if an additional concentration ratio is used at reception (in addition to the 4:1 ratio specified by NEMA for P5), <b>this additional concentration ratio can be used for P6 and KP5.</b>	<b>Minimal sphere size</b> has to be the same (at least) as the one observed during reception test P6.
KP6	Image homogeneity and precision of quantification	Half-yearly	See Annex 6 for detailed instructions.  <b>The coefficient of variation and mean background activity are measured and compared with given tolerance intervals.</b>	All	Radionuclide: F-18 with 50-100 MBq. Use of a cylindrical fillable phantom (same as for P7).  For the centers undergoing accreditation tests for quantitative measurements (e.g., EARL), this test can be skipped, and the accreditation test can be presented.	Tolerance intervals are given in Annex 6.

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




Figures taken from:  
Jonathan I Gear *et al.* Radioactive 3D printing for the production of molecular imaging phantoms. 2020 Phys. Med. Biol. 65 175019 DOI 10.1088/1361-6560/aba40e

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

Printing with [ $^{18}\text{F}$ ]-FDG from aqueous solution on a masked stereolithography (MSLA) 3D printer (SL1, Prusa Research, Prague, Czech Republic).

Figures and text taken from:  
 Gillett, D., Marsden, D., Ballout, S. et al. 3D printing 18F radioactive phantoms for PET imaging. EJNMMI Phys 8, 38 (2021).  
<https://doi.org/10.1186/s40658-021-00383-6>

## Radioactivity leaching

“The amount of radioactivity that leached into the water of the phantom was **0.72%** of the activity in the helix. This was calculated by taking a sample from the water at 3 h after the phantom (Fig. 6d) being in the water.”

# Results: $^{68}\text{Ge}$ -extraction, SLA 3D printing and source tightness

## 1. $^{68}\text{Ge}$ -extraction

Radioactivity conc. ratio  
organic/water: >3700:1

## 2. Stereo-lithography


3D print,  
 $\text{Ge-68}$

## 3. Safety: Wipe Test


- 5x Water: Below LOQ
- 5x Ethanol: Below LOQ

## 4. Enhanced Safety Test: Washout at 50°C, 1h


# Results: PET/CT and homogeneity measurements

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
Line phantom for measuring spatial resolution of clinical SPECT systems. Glass capillaries are currently filled with 55 to 80 MBq Technetium-99m solution (a) and imaged on the SPECT system laid out in a square pattern (b). Line spread function can then be derived from the images for two spatial directions and for different locations in the field of view. (IAEA Human health series No. 6, Quality assurance for SPECT systems, IAEA, 2009). We are planning to replicate this phantom with long lived Cobalt -57 using our manufacturing methods. Mouse spine phantom print model (c) and mouse head phantom (d) with **S**kull, **B**rain, and **T**umor for preclinical molecular imaging. These types of phantom will be either printed with Germanium-68 or Cobalt-57 containing resin.

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